

in vitro Transcription**REAGENTS**

- T3 RNA polymerase (P2083, Promega)
- UltraPure DNase/RNase-Free Distilled Water (10977, GIBCO)
- 10 mM NTP (Invitrogen)
ATP: 18330-019, GTP: 18332-015, CTP: 18331-017, UTP: 18333-013
- m⁷G(5')ppp(5')G RNA Capping Analog (15619-018, Invitrogen)
- Recombinant RNase Inhibitor (2311A, TaKaRa)
- TURBO DNase (AM2238, Ambion)

REAGENT SETUP

- 1) Dilution of Capping Analog to 10 mM by UltraPure Water
- 2) Preparation of NTP mix

NTP mix for CAS9 mRNA synthesis (NTP with CAP analog)

ATP:CTP:UTP:GTP:GTP capping analog = 4:4:4:1:3 (molar ratio)

NTP mix for guideRNA synthesis (NTP without CAP analog)

ATP:CTP:UTP:GTP = 1:1:1:1 (molar ratio)

PROCEDURE

- 1) Incubate following mix at 37°C for 2 hr.

	(μl)
5x buffer	8
NTP mix	6.4
0.1 M DTT	4
Linearized Template DNA	3000 ng of CAS9 template, or 1000 ng of gRNA template
RNase inhibitor	2
Polymerase	4
UltraPureWater	up to 40
<hr/>	
Total	40 μl

- 2) Add 2 μl of DNase, and mix gently by pipetting.
- 3) Incubate at 37°C for 30 min.
- 4) Add 60 μl of UltraPure Water and 100 μl RNase-free PCI, and mix by vortexing.
- 5) Centrifuge at 20000 x g for 5 min, and transfer the supernatant to a fresh tube.
- 6) Add 40 μl of 10 M of ammonium acetate and 370 μl of ethanol. Mix by vortexing.
- 7) Keep at -80°C for 15 min.
- 8) Centrifuge at 20000 x g for 30 min. Discard sup and add 200 μl of 70% EtOH.
- 9) Centrifuge at 20000 x g for 5 min. Discard sup.
- 10) Dilute RNA to 200 μg/ml by 10-30 μl of UltraPure Water. Aliquot and store at -80°C.